

Chromatography Vial



Product Categories: [Chromatography Vials Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/chromatography-vial/>

Product Description

Chromatography Vial - Laboratory Glassware Equipment

Designed to use separating organic and inorganic compounds for analysis, chromatography vials come in different finishes to guarantee chemical compatibility. Auto sampler compatible closures come in different colors, sizes, and septa materials. Offered with a write-on patch, possible identification errors are eliminated. Supporting use with inserts or already possessing a conical interior shape, maximum retrieval of contents to eliminate waste. Select from white, amber, or clear for trouble-free mixture and solvent observation.

Features:

Long lasting performance

Durability

Lightweight

Accurate dimensions

Sleek design

Robust construction

Very high durability

Longer functionality

Auto-zero on Reservoir Plain with Stop-cock

Screw Cap, made from clear glass tubing with Silicon 1 mm Septa & Cap.

Capacity Approx -4 ml, 6 ml, 10 ml, 20 ml.