


Clinical Thermometer



Product Categories: [Chemistry Equipments](#), [Thermometer & Hygrometer](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/clinical-thermometer/>

Product Description

Clinical Thermometer

Clinical thermometer designed to measure the temperature of the human body. Graduated to cover a range a few degrees on either side of the normal body temperature.

Prismatic Capillary, Chemical thermometer.

We hold expertise in exporting and importing Thermometer Clinical, which are made using highest quality both in material and design.

These are absolutely reliable, safe and accurate to use, clinical thermometer comes in various measuring ranges. Moreover, the offered thermometer is made with permanent fused scale, with constriction in bore to prevent mercury from going back.

Feature

Prismatic Capillary

Made from japan glass

With perfect pressure detection and measuring accuracy