

Co-Axial Binocular Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/co-axial-binocular-microscope/>

Product Description

Co-Axial Binocular Microscope - Scientific Equipment

Coaxial Binocular Microscope comprising of large stand and binohead, coaxial mechanism for coarse and fine motions, coaxial mechanical stage, quadruple nose-piece with ball-bearings, built-in solid state transformer for 6 V/20 watt Halogen lamp binocular head fitted with superior hard coated stain free prisms. The product is nicely packed in wooden box.

Nose piece: Reverse type Quadruple ball bearing nose-piece. Binocular head 360* rotate able and 45* inclined head with all anti-fungus coated prisms.

Focusing system: Co-axial focusing system with three planetary gears assb. With slide protection lock.

Eyepieces & objectives: Set of 4 x, 10 x, 40 x, and 100 x and a paired set of 10 x wide field eyepieces.

Illumination: 6 v, 20 w halogen lamp needs 220 v (50 Hz), 110 V(60 Hz) optional supply.

Optional: LED illumination system also available with battery backup.

Specifications:

Objectives: 4 x, 10 x, 40 x, 100 x oil immersion

Eyepieces: 10 x (WF)