

## Comparative Germination Model



**Product Categories:** [Educational equipments](#), [Educational Models](#), [Zoology Models & Botany Model](#)

**Product Tags:** [biology lab](#), [botany model](#), [fibre model](#), [science model](#), [zoology](#)

**Product Page:**

<https://www.labappara.com/product/comparative-germination-model/>

### Product Description

#### **Comparative Germination Model- Scientific Educational Model**

Dimensional three different types of plants on one model. Compare monocot, dicot, and gymnosperm germination with this one high-quality, exceptionally detailed model. The monocot type is represented by a rye shoot enlarged 10x. The dicot by a two-leaf bean shoot enlarged 5x, and the gymnosperm by a star-shaped fir leaflet enlarged 20x. Students can be study root development, the methods by which the different plants sprout, and the structure of each type of seedling because each component of the model is removable for closer examination. The eight-piece model is mounted on a base and comes with a study guide identifying 9 structures.