

COMPOUND METAL STRIP (EXPORT)



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/compound-metal-strip-export/>

Product Description

Compound Metal Strip - Physics Equipment

Compound Bars are also known as Bimetallic Strips.

For the experiments on thermal expansions in physics labs.

Multiple metal combinations turn a simple activity into an advanced learning experience.

Two Different Metals riveted together end to end with wooden handle.

Compare the behavior of four combinations.

Brass & Aluminum

Iron & Brass

Iron & Copper

Brass & Copper

make conclusions about appropriate applications of each metal.