

# COMPOUND METAL STRIP



**Product Categories:** [Heat And Thermodynamics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/compound-metal-strip/>

## Product Description

Compound Metal Strip - Physics Equipment

Compound bars are also known as Bimetallic Strips.

For the experiments on thermal expansions in physics labs.

Multiple metal combinations turn a simple activity into an advanced learning experience.

Two Different Metals riveted together end to end with wooden handle.

Compare the behavior of four combinations: Brass/Aluminum, Brass/Iron, Copper/Iron, and Brass/Copper, and ask students to make conclusions about appropriate applications of each metal.

Brass & Aluminum

Iron & Brass

Iron & Copper

Brass & Copper