

Concave Mirror



Product Categories: [Lenses](#), [Optical Instruments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/concave-mirror-2/>

Product Description

Concave Mirror - Scientific Laboratory Educational Equipment

A concave mirror, or converging mirror, has a reflecting surface that bulges inward (away from the incident light). Concave mirrors reflect light inward to one focal point. They are used to focus light. Mirrors, highly slivered, optically true (spray painted at the back) with a beautiful ring on its outer surface

Diameter	Focal length
-----------------	---------------------

50mm	10-100cm
------	----------

75mm	10-100cm
------	----------