

Condensor Universal Clamp



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/condensor-universal-clamp/>

Product Description

Condensor Universal Clamp - Scientific Equipment

Condenser may refer to a device that is used to condense a fluid from its gaseous state to liquid.

Feature -

Brass oxidised with boss head.

Made of special aluminium Alloy.

Provide with 7" mild steel rod C.P.

And holds upto 75 mm dia articles.

Clamp is staple painted, jaws are cork lined.

Jaws are made from Aluminium casted material.

Jaws are covered with thick rubber sheet.

Able to catch up to 75 mm dia objects.

150/200/250mm long chrome plated rod.

Duly powder coated finishing.