

## Cone Borosilicate Glass Imhoff Cone



**Product Categories:** [Cones Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:**

<https://www.labappara.com/product/cone-borosilicate-glass-imhoff-cone/>

### Product Description

Cone Borosilicate Glass Imhoff Cone - Laboratory Glassware Equipment

Glass Construction for determination of small amount of Suspended matter in Sewage by the Imhoff sediment Method. Cones, Imhoff, Sediment, Shape Tip:

Graduated from 0 to 1 ml in 0.1 ml division

Graduated from 1 to 10 ml in 0.5 ml division

Graduated from 10 to 40 ml in 1 ml division

Also marked at 1000 ml.

Approx Total Height 454 mm

Approx Top Dia 100 mm Cones, Imhoff, Sediment, Blunt Tip:

Graduated from 0 to 1 ml in 0.1 ml division

Graduated from 1 to 10 ml in 0.5 ml division

Graduated from 10 to 40 ml in 1 ml division

Also marked at 1000 ml.

Approx Total Height 439 mm

Approx Top Dia 100 mm

Measure the volume of settled in water and sewage

Shatter proof poly carbonate cone is glass clear.

Graduations up to 1000 ml for samples with a high volume of settled matter.

Waterlight polyethylene threaded allows easy withdrawal of simplifies cleaning.

**Uses:**

Used to measure the volume of settle able solids in a specific volume.

**Features:**

Scratch-proof

Excellent finish

Durability