

## Cones, Imhoff, Sediment, Blunt Tip



**Product Categories:** [Cones Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:**

<https://www.labappara.com/product/cones-imhoff-sediment-blunt-tip/>

### Product Description

Cones, Imhoff, Sediment, Blunt Tip - Laboratory Glassware Equipment

1 liter blunt tip cones are for use with heavy sediments, where ability to read small volumes is unimportant. The blunt tip increases its ruggedness and cleaning ease.

Height: approximately 435 mm. Top Inner Diameter: approximately 110 mm.

Graduation intervals: 0 ml to 1 ml in 0.1 ml divisions; 1 ml to 10 ml in 0.5 ml divisions; 10 ml to 40 ml in 1 ml divisions. Also marked at 1000 ml. Material: Borosilicate Glass

Cones, Imhoff, Sediment, Blunt Tip Graduated from 0 to 1 ml in 0.1 ml.

Division, 1 to 10 ml. in 0.5 ml. and 10 to 40 ml. in 1 ml.,

Approx. Top Dia 100 mm.