


## Cones, Imhoff, Sediment, Sharp Tip



**Product Categories:** [Cones Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:**

<https://www.labappara.com/product/cones-imhoff-sediment-sharp-tip/>

### Product Description

Cones, Imhoff, Sediment, Sharp Tip - Laboratory Glassware Equipment

Cones, Imhoff, Sediment, Sharp Tip Graduated from 0 to 1 ml in 0.1 ml.

Division, 1 to 10 ml. in 0.5 ml. and 10 to 40 ml. in 1 ml.

Also marked at 1000 ml. Approx.

Total Height 454 mm.,

Approx. Top Dia 100 mm.

### Features:

Durable

Light weight

Flawless finish