

## Conical Flask



**Product Categories:** [Laboratory Plastic Ware](#), [Liquid-Measurement](#)

**Product Tags:** [Laboratory Equipment](#), [Laboratory plastic ware](#)

**Product Page:** <https://www.labappara.com/product/conical-flask-2/>

### Product Description

Conical Flask - Scientific Laboratory Equipment

Conical Flasks are made up of polypropylene, these cone shaped flasks are rigid, translucent & autoclavable. The screw cap keeps the contents of the Conical Flasks contamination free and makes the flask leak proof. The threads of the Conical Flasks should be disengaged before autoclaving.

#### **Features:**

Superior strength

Perfect finish

Leak proof

Resistant to abrasion

**capacity :** 100, 250, 500 ml