

## Conical Flask



**Product Categories:** [Bottles And Assesories](#), [Laboratory Plastic Ware](#)

**Product Tags:** [Labaoratory Equipments](#), [Laboratory plastic ware](#)

**Product Page:** <https://www.labappara.com/product/conical-flask/>

### Product Description

Conical Flask - Scientific Laboratory Equipment **Conical Flasks** are made up of polypropylene, these cone shaped **flasks** are rigid, translucent & autoclavable. The screw cap keeps the contents of the **Conical Flasks** contamination free and makes the **flask** leak proof. The threads of the **Conical Flasks** should be disengaged before autoclaving.