

Conversion Of Galvanometer Into Ammeter



Product Categories: [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

Product Page:

<https://www.labappara.com/product/conversion-of-galvanometer-into-ammeter/>

Product Description

Conversion Of Galvanometer Into Ammeter

Features :

Instrument comprises of DC Regulated Power Supply 0-5VDC/150mA, two round meters (One milliammeter & One Galvanometer), Three set of Shunts resistance mounted behind the panel, connections of supplies, Meters & Shunts brought out at 4mm Sockets.