

## Cook's Pole Climbing



**Product Categories:** [Pharmaceutical Equipments](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/cooks-pole-climbing/>

### Product Description

#### Cook's Pole Climbing - Pharmaceutical Equipment

This is an improved user friendly model with both light & sound stimuli can be selected by a switch. Shock stimulus is of 200 volt 2 M amp 50 HZ A.C. Application of shock can be either manual or automatic and also ready or delayed. Both duration of delay and application of shock are controlled by 11 position timer switch.

Cook's Pole Climbing Apparatus use to study cognitive function, mainly a response to conditioned stimuli during learning & its retention.

The apparatus has an experimental with the floor grid in a soundproof enclosure.