

# Coplin Jar



**Product Categories:** [Laboratory Plastic Ware](#), [Slide Storage](#)

**Product Tags:** [biology lab equipment](#), [Laboratory Equipments](#), [Laboratory plastic ware](#)

**Product Page:** <https://www.labappara.com/product/coplin-jar/>

## Product Description

Coplin Jar - Scientific Laboratory Plastic Ware

Coplin Jars are made up of polypropylene and can hold 10 sliders of the size 1" x 3" back to back. The interior of the Coplin Jar is grooved in way to hold slides in erect position. Domed & shallow threaded screw cap provides protection to the slides.

**Uses:**

Used for holding microscope slides during staining procedures.

**Features:**

Superior strength

Impeccable finish

Highly durable

Particulars	Packing
Coplin Jars	12 Pc