

Coplin Jar



Product Categories: [Laboratory Plastic Ware](#), [Slide Storage](#)

Product Tags: [biology lab equipment](#), [Laboratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/coplin-jar/>

Product Description

Coplin Jar - Scientific Laboratory Plastic Ware

Coplin Jars are made up of polypropylene and can hold 10 slides of the size 1" x 3" back to back. The interior of the Coplin Jar is grooved in way to hold slides in erect position. Domed & shallow threaded screw cap provides protection to the slides.

Uses:

Used for holding microscope slides during staining procedures.

Features:

Superior strength

Impeccable finish

Highly durable

Particulars

Packing

Coplin Jars

12 Pc