

## COPPER CALORIMETER COMPLETE



**Product Categories:** [Heat And Thermodynamics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/copper-calorimeter-complete/>

### Product Description

Copper Calorimeter Complete - Physics Instrument

Used for measuring the heat of chemical reactions or physical changes as well as heat capacity.

It consists of a Copper vessel housed in a rectangular wooden box fitted with a thermometer holder

Vessel of Copper in outer wooden case

With thermometer holder and stirrer, well insulated in card box.

Used to measure specific heats of substances

Its simple construction makes it easy to use while delivering accurate results.

The space between the Copper vessel and the housing is packed with felt for thermal insulation.

Supplied with a Copper stirrer, but without thermometer.

Available in seasoned wood ,and Teak wood

Supplied according to the demand of customer

Seasoned Wood 4" x 3"

Seasoned Wood 3" x 2"

Teak Wood 3" x 2"

Teak Wood 4" x 3"