


Copper Colorimeter (Teak Wood)



Product Categories: [Analytical Lab Instruments](#), [Calorimeters](#)

Product Tags: [science lab equipment](#)

Product Page: <https://www.labappara.com/product/copper-colorimeter-teak-wood/>

Product Description

Copper Colorimeter (Teak Wood) - Analytical Laboratory Equipment

Vessel of copper in outer wooden case with thermometer holder and stirrer, well insulated in card box.

It consists of a rectangular teak wood box of fitted with a thermometer holder and packed with felt for thermal insulation.

Supplied with a copper stirrer but without thermometer.

This calorimeter is domestic kcal instrument used to measure heat energy. It consists of a rectangular teak wood box of size 9.5"cm fitted with a thermometer holder and packed with felt for thermal insulation.

Features:

Sturdy construction

Longer service life

Superior performance

It consists of a rectangular teak wood box of size 3"x 2" , 4" x 3" fitted with a thermometer holder and packed with felt for thermal insulation.

Supplied with a copper stirrer but without thermometer.

Available in size - 3"x 2" , 4" x 3"