

## Copper Colorimeter



**Product Categories:** [Analytical Lab Instruments](#), [Calorimeters](#)

**Product Tags:** [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/copper-colorimeter/>

### Product Description

Copper Colorimeter - Analytical Laboratory Equipment

Vessel of copper in outer wooden case with thermometer holder and stirrer, well insulated in card box.

It consists of a rectangular wood box.

Fitted with a thermometer holder and packed with felt for thermal insulation.

Supplied with a copper stirrer but without thermometer.

- Available in size - 3"x 2" , 4" x 3"