

Copper Colorimeter



Product Categories: [Analytical Lab Instruments](#), [Calorimeters](#)

Product Tags: [science lab equipment](#)

Product Page: <https://www.labappara.com/product/copper-colorimeter/>

Product Description

Copper Colorimeter - Analytical Laboratory Equipment

Vessel of copper in outer wooden case with thermometer holder and stirrer, well insulated in card box.

It consists of a rectangular wood box.

Fitted with a thermometer holder and packed with felt for thermal insulation.

Supplied with a copper stirrer but without thermometer.

- Available in size - 3" x 2" , 4" x 3"