

# Copper Turning



**Product Categories:** [Chemistry Equipments](#), [General Chemistry Labware](#)

**Product Tags:** [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/copper-turning/>

## Product Description

Copper Turning - Chemistry Equipment

Copper Turnings are the shavings of copper metal (often left after manufacturing of copper goods in factories are used mostly in schools for studying chemical reactions as it provide more surface area for reaction than a piece of copper ).

They have the same color as copper (reddish brown).

Features:

Reliable usage

Recyclable

Purity