

Copper Voltameter



Product Categories: [Physiology Apparatus](#)

Product Tags: [Labaoratory Equipment](#)

Product Page: <https://www.labappara.com/product/copper-voltameter/>

Product Description

Copper Voltameter - Laboratory Equipment

For experiments in electrolysis 2 thick copper plates of 70 x 50 mm and a central plate 70 x 45 mm forming the cathode are suspended from bakelite cover in a glass jar 150 x 100 mm (height x dia). The outer plate from anode and the central one is cathode. The central plate removable.

The S.I of quantity of electricity is the coulomb.

The Voltmeter should not be confused with a voltmeter, which measures electric potential.

The SI unit for electric potential is the volt.

Copper voltameter and with jar.