

Cork Borer



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/cork-borer/>

Product Description

Cork Borer - Laboratory Equipment

A cork borer, often used in a chemistry or biology laboratory,

A metal tool for cutting a hole in a cork or rubber stopper to insert glass tubing.

Cork borers usually come in a set of nested sizes along with a solid pin for pushing the removed cork (or rubber) out of the borer.

Made of plated sharpened brass tubing with attached handles including a steel punch.

Feature -

Made from thick brass.

Set of 6 to 12 borers in velvet box.