

# Crompton Potentiometer



**Product Categories:** [Electrical Bridge's](#), [Electrical Bridge's](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/crompton-potentiometer-2/>

## Product Description

Crompton Potentiometer

**Features:**

Rust proof body

Durability

Easy to operate

**Specification:**

Direct reading type, improved from with standardizing coil consisting of 8 coils and a circular slide wire

The scale is divided into 200 parts with potential drop of 200 mill volts

Fine and coarse rheostats are provided for adjusting the slide wire current

Range 0.001 to 1.8 volts, managing coils, provided with or without standard cell, Galvanometer and Rheostat