

Cross Connectors



Product Categories: [Connectors And Stirrers](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/cross-connectors/>

Product Description

Cross Connectors - Scientific Laboratory Equipment

These four way connectors proved to be an excellent option for achieving leak proof joints in tubing connections. The connectors are basically pipe-like structures in shapes with ringed edges. These are used for connecting tubes to each other.

The ringed edged prevent leakage & slippage

Material: Polypropylene

Particulars	Packing
Cross Connector - 6 mm	36 Pcs
Cross Connector - 8 mm	36 Pcs
Cross Connector - 10 mm	36 Pcs