


Crucible Tongs C.P



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/crucible-tongs-c-p/>

Product Description

Crucible Tongs C.P - Chemistry Laboratory Equipment

It is available in different sizes & specifications at competitive prices.

Material:-

High Grade Raw Material

Cutting Edge Technology

Feature:-

Crucible tongs are scissor-like tools, but without blades.

These tools are replaced with two pincers or pieces of metals that concave together,

Which allow the users to grasp a hot crucible, flasks, evaporating dishes, or even small beakers .

They are made of durable metals, allowed them to endure a very hot temperature when touching with the hot objects.

Sleek Design

Comfort Grip

Rust Resistant

Available in size - 6", 8" (Iron & S.S)