

## Crucible



**Product Categories:** [Laboratory Glassware](#), [Sintered Ware Borosilicate Glass](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/crucible/>

### Product Description

Crucible - Laboratory Glassware Equipment

Crucible gooch is a filtration device for laboratory use. It can also be used for collection and processing of biological tissue samples within the same container. Our crucible gooch is widely used in various labs, institutes and colleges. Crucible Gooch type with Sintered Disc of Porosity G-1 or G-2 or G-3 or G-4.

Features:

Durability

Longevity

Quality assured

Cat. No.	Capacity
328/1	15
328/2	30
328/3	50