

CRYO Box (PP)



Product Categories: [Laboratory Plastic Ware](#), [Vials](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/cryo-box-pp/>

Product Description

CRYO Box (PP) - Scientific Laboratory Equipment

Moulded in polypropylene, the top of the box is attached to the bottom through well designed built-in hinges. The bottom has built - in spaces to hold vials up to 2 ml. capacity. Available in two different sizes, these autoclavable boxes bear alpha numeric index for easy identification of different vials placed inside.

Material : Polycarbonate

Particulars

Packing

81 Places for 1 ml. or 1.8 ml Cryo Vials

4 Pcs

100 Places for 1 ml or 1.8 ml Cryo Vials

4 Pcs