

## CRYO Rack



**Product Categories:** [Laboratory Plastic Ware](#), [Vails](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/cryo-rack/>

### Product Description

CRYO Rack - Scientific Laboratory Equipment

Cryo Rack is moulded in Polycarbonate making it tough, durable & autoclavable. This handy rack can hold 50 vials at a time. The unique feature of this rack is that vials can be unscrewed without actually picking them up as the vial will lock in every well and hence will not twist. This rack is provided with four anti skid rubbers & each well is identified with an alpha-numeric index.

#### Particulars

50 places for 1, 1.8 & 4.5 ml Cryo Vials

#### Packing

2 Pcs