

## CRYO Vial-Internal Thread



**Product Categories:** [Laboratory Plastic Ware](#), [Vials](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/cryo-vial-internal-thread/>

### Product Description

CRYO Vial-Internal Thread - Scientific Laboratory Equipment

Adding to the existing range of Cryo Vials which has internal thread in it. The cap is fitted with an O-ring screws in to the vial and makes it a leakproof Cryo Vial. These Vials are provided with white, prominent graduation and a long marking area for sample identification.

**Material :** Polypropylene

<b>Particulars</b>	<b>Packing</b>
--------------------	----------------

2.0 ml	500 Pcs
--------	---------