

CRYO Vial



Product Categories: [Laboratory Plastic Ware](#), [Vials](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/cryo-vial/>

Product Description

CRYO Vial - Scientific Laboratory Equipment

Cryovial are normally used for storage of Biological material, human & animal cells. These sterilized tubes are provided with white prominent graduation and a long marking area for sample identification. The cap of these Cryovial is provided with a seal ring inside to make it completely leak proof. The longitudinal grooves in the base give a non-twisting feature to these vials when placed in a cryo rack.

Particulars	Packing
1.0 ml	1000 Pcs
1.8 ml	1000 Pcs
4.5 ml	500 Pcs