

CUBES OF METAL



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/cubes-of-metal/>

Product Description

Cubes Of Metal - Physics Equipment

Set of 6 different metal cubes

These metal blocks are intended for use as sinkers or for specific gravity measurements.

Set of six different metals. The set consists one piece each of brass, aluminum, copper, zinc, iron and lead.

These are used for measurement of density and specific gravity.

Features:

Sturdiness

Fine finish

Perfect strength

SIZE:-

10 mm

20 mm

25 mm