

CYLINDER OF METALS



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/cylinder-of-metals/>

Product Description

Cylinder Of Metals - Physics Equipment

Set of 6 different metal cylinders

These cylinders are often used in exercise with Vernier and micrometer calipers as well as for density and specific gravity measurement.

The ends are carefully turned and surfaces are smooth.

Set of six different elements in cylindrical form for density in electrical and thermal experiments.

The set consists of one piece each of Brass, aluminum, copper, zinc, iron and lead solid cylinders.

Features:

Sturdiness

Fine finish

Perfect strength

SIZE:-

25 x 12 mm

37 x 12 mm

50x12 mm