

Cylindrical Magnet



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/cylindrical-magnet/>

Product Description

Neodymium Cylindrical Magnet - Physics Instrument

Cylindrical Magnet is designed using optimum grade raw material and superior technology with the help of our highly skilled team.

Made up of 12 mm dia chromium steel bars red and blue painted with keepers.

Magnetizing and handling:

Neodymium magnets are very brittle and possess very strong magnetic field.

These magnets can chip or even can cause injuries if not properly handled.

It is suggested to magnetize these magnet after assembling in machines

Applications:

Laptop computer disk drives, printers, speakers, novelties, toys, DC motors and magnetic separators

Properties:

It possesses high energy for small sizes.

It has low working temperature.

Features:

Strong construction

Extremely strong

Excellent magnetic property

Precise dimensions

High efficiency

Good Quality

Available in size -25 mm, 50 mm, 75 mm, 100 mm, 125 mm, 150 mm