

## Dale Bath Apparatus



**Product Categories:** [Pharmaceutical Equipments](#)

**Product Tags:** [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/dale-bath-apparatus/>

### Product Description

#### **Dale Bath Apparatus - Pharmaceutical Equipment**

A Perspex organ bath with fine thermostatic temperature control, electronic stirrer and complete set of accessories, making it a versatile unit for all research work. Bath 9"×6"×8", with leak proof dovetailed joints, resting in sturdy metal frame; top of the bath carries a 1"×3/8" hick band, for fixing accessories, inside the bath is a metal sheathed heating element (easily replicable) controlled by means of a fine adjustable thermostat with sensitivity  $\pm 1^{\circ}\text{C}$ . The stirrer unit is completely independent and dips in the bath in one corner, fixing on its special stainless steel upright. A M.S. Junction box is fitted to the bath frame and carries all the electrical leads and connections two pilot lamps indicate functioning of the thermostat and stirrer.

Standard components are:

Frontal Writing Lever-1

Organ Tube-1

Aeration Tube-1

Stirrer-1

Stainless Steel capillary levers

The organ bath is a traditional experimental set-up that is commonly used to investigate the physiology and pharmacology of in vitro tissue preparations.