

Davies Condensers



Product Categories: [Condensers Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/davies-condensers/>

Product Description

Davies Condensers - Laboratory Glassware Equipment

Highly efficient condenser with large cooling surface.

Inner cooling surface gives 'baffle' effect.

Outer cooling surface prevents 'creep' with low boiling point liquids.

With detachable plastic screw thread connectors for easy and safe attachment of flexible tubing.

Key Features:

Impact resistant

Resistant to high heat

Longer life

Available in effective length (mm) - 150 mm, 200 mm, 250 mm, 300 mm, 400 mm, 500 mm, 600 mm.