

Deflagrating Spoon (Iron)



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/deflagrating-spoon-iron/>

Product Description

Deflagrating Spoon (Iron) - Scientific Equipment

Deflagrating-spoon (chemistry) A long-handled spoon with a cover. Lowered into a glass vessel filled with a gas to demonstrate deflagration.

A deflagration spoon is used in chemistry experiments when the burning of [phosphorus](#), [sulfur](#) or another substance is required.

It is capable for high temperatures.

It is equipped with features to facilitate working in this environment.

Such as a long handle angled to keep heat away from the hand.

AVAILABLE IN:

Iron wire brass disc.

Export quality