

## Deflection Compass (Permanent Magnet)



**Product Categories:** [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/deflection-compass-permanent-magnet/>

### Product Description

#### Deflection Compass (Permanent Magnet) - Physics Equipment

Compressing as thick Bakelite moulded case of 125 mm dia. with black aluminum dial 0 graduated 0 - 90° four times, anti parallax mirror slot.

The special permanent magnet (small in size) and a light maximum damping) supported on a steel pivot by jewel bearing sets.

The permanent magnet is designed in such a way so that it has maximum magnetic fields strength.

Compass Deflection for Physics Lab - We are one of the leading exporters who offers a full line of Laboratory Physics Products.

Physics Products as per technical specification and tender requirements for labs, educational institutions & universities.

#### **Specifications:**

Accurate result

Available as per clients' demands

Easy to use

Precision-designed

Low maintenance

Durability