

Desauty Bridge with Head Phane & Audio Oscillator 1 KHz



Product Categories: [Electrical Bridge's](#), [Electrical Bridge's](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/desauty-bridge-with-head-phane-audio-oscillator-1-khz/>

Product Description

Desauty Bridge with Head Phane & Audio Oscillator 1 KHz

Features:

- Compact design
- High performance
- Superior function

De-Sauty Bridge for capacitance measurement:

- Complete in all respect, no external component required
- De-Sauty Bridge to determine Capacitance of a Capacitor
- R1 3 decade resistance dials having range X 10, 100 & 1K ohms
- R2 Three more decade of same value as above
- C1 Standard Capacitors having .1 mfd and .2 mfd
- C2 Four unknown Capacitors are provided

Specifications:

Inbuilt A.C supply of frequency 1 KHz and Digital Null Detector for AC balance is provided. 6 interconnection Leads with 4mm banana plugs, 4 unknown Capacitances are provided on board. Comprehensive User's Manual will explain required results.