

## Desiccator Plain



**Product Categories:** [Desiccators Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/desiccator-plain/>

## Product Description

Desiccator Plain - Laboratory Glassware Equipment

Desiccator Plain with knobbed lid, Ground Flanges Neutral Glass, with Perforated Porcelain Plate.

Plain Desiccators are physically the same as vacuum desiccator. Only difference is the absence of stop cock. Grease should be applied on the flanges to make it air tight.

Available in diameter(mm) - 160mm, 210mm, 240mm, 300mm, 120mm