

## Desiccator, Vacuum



**Product Categories:** [Desiccators Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/desiccator-vacuum/>

### Product Description

Desiccator, Vacuum - Laboratory Glassware Equipment

Desiccator, Vacuum With tubular Lid, Ground Flanges with Glass stopcock made of Neutral glass, supplied with Perforated Porcelain Plate.

Desiccators are seal-able enclosures containing desiccants used for preserving moisture-sensitive items such as cobalt chloride paper for another use. A common use for desiccators is to protect chemicals which are hygroscopic or which react with water from humidity.

#### **Features:**

Long working life

Clear body

Not compatible with organic solvents or vapors.

Room temperature use only.

Precisely designed

Optimum functionality

Available in dia(mm) - 160 mm, 210 mm, 240 mm, 300 mm