

Digital Conductivity Meter



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Equipment](#)

Product Page: <https://www.labappara.com/product/digital-conductivity-meter-2/>

Product Description

Digital Conductivity Meter - Analytical Equipment

Three and half digit LED display with facility for temperature compensation and digital cell.

constant adjustment along with conductivity cell & other accessories.

Highly Stable and Accurate Digital Conductivity Meter

Cell Constant Adjustment Facility is provided in this Digital Conductivity Meter

This Digital Conductivity Meter available in Lab & Field models

Temperature Coefficient Correction facility also provided in the Digital Conductivity Meter

Range

0 - 1999 pS/m

Resolution

1 pS/m

Accuracy

2% of reading \pm 1 pS/m

Over-RangeIndicator

"1" on left side of display

Calibration

Internal source, field adj.

Controls

2 Pushbuttons, "M" (measure)& "C" (calibration)

Display

Liquid crystal, 3 ½ digits, ½ inch (.27 cm) high

SafetyGround

Banana jack on electronics assy.

Power

Battery, 3 each, 6 volt alkalineOp Temp

Range

32° to 165° F (0° to 75° C)

CarryingCase

Hard, solvent resistant,molded plastic