

Digital Dissolved Oxygen Meter



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Lab Equipment](#)

Product Page:

<https://www.labappara.com/product/digital-dissolved-oxygen-meter-2/>

Product Description

Digital Dissolved Oxygen Meter- Analytical Lab Equipment

Three and half digit LED display with gold/silver probe.

LCD display ABS plastic body, along with Gold/Silver DO Probe.

9V battery operated.

a) Dissolved Oxygen Probe (Amperometric Gold/Silver membrane Type)

b) A set of 10 DO membranes

Specifications:

Range 0 to 20.0 ppm

Resolution 0.1 ppm

Accuracy ± 0.2 ppm

Power Requirement 220V AC

Complete with Instruction Manual & Dust Cover.