

## Digital Dissolved Oxygen Meter



**Product Categories:** [Lab Measuring Instruments](#)

**Product Page:**

<https://www.labappara.com/product/digital-dissolved-oxygen-meter/>

**Product Description**

### **Digital Dissolved Oxygen Meter - Scientific Equipment**

3 and half digit LED display with gold/silver probe.

Range 0 to 20.0 ppm

Resolution 0.1 ppm

Accuracy  $\pm 0.2$  ppm

Power Requirement 220V AC

Complete with Instruction Manual & Dust Cover.