

Digital Friability Test Apparatus (Double Drum)



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/digital-friability-test-apparatus-double-drum/>

Product Description

Digital Friability Test Apparatus (Double Drum) - Analytical Laboratory Equipment
Measures obsession strength of tablets with programmable timer.

Rotation speed 25 rpm, accuracy ± 1 , auto & manual mode of operation with
Double acrylic drums.

Results are display and it has a 3 digit presettable control.

- a) Acrylic Tray for Single Drum
- b) Acrylic Tray for Double Drum
- c) Drum locking Knob

Model

Instrument

Specification

LT-70

Digital Friability Test Apparatus (DOUBLE DRUM)

Measures Obsession Strength of Tablets with programmable timer. Rotation Speed 25 rpm +/- 1, auto & manual mode of operation with 2 acrylic drums, Results are displayed on 3 digit LED display and it has 3 digit prestable control.