

DIGITAL MELTING POINT APPARATUS



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Lab Measuring Instrument](#), [Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/digital-melting-point-apparatus-4/>

Product Description

Digital Melting Point Apparatus - Analytical Laboratory Instrument
Silicon Oil Bath, Built-in Magnetic Stirrer with electronic speed controller.
Electronic Controller for adjusting the heating rate.
Digital Display of Temperature with PT-100 Sensor.
Glare-free back ground light with adjustable light Intensity.
Provision for holding the melting point display Suitable to determine boiling point also and accuracy 1.0 Degree.

Features:

- Precisely designed
- Long working life
- Low power consumption

Specifications

Display

Digital 3 1/2 Digit, Bright LED

Accuracy

1%

Operation

Fully Automatic

Dimensions

L 405 x B 360 x H 235 mm (approx)

Weight

5 Kg (approx.)

Power

230V + 10%, 50Hz

Accessories

Capillary Tube - 100 No.

Cleaning Tool - 1 No.