

Digital PH Meter



Product Categories: [Laboratory Equipments](#), [Ph Meters](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/digital-ph-meter/>

Product Description

Digital PH Meter - Scientific Laboratory Equipment

PH meters are stable test instruments used to measure acidity or alkalinity of a solution. Finding extensive usage in chemical, medical, biology and physics laboratories these meters have been provided with advanced features for enhanced performance. Our pH meters have digital LED screen which displays the exact pH value of a solution promptly.

Features:

Robust design

Optimum quality

Easy to install

Flawless performance

Easy usage

Long service life

Specifications:

Display 3½ LED

Range 0-14 pH 0±1999mV

Resolution 0.01 pH 1mV

Temperature Compensation Manual Automatic

Power Requirement 220V AC

With Dust Cover & Instruction Manual.