

## DIGITAL POTENTIOMETER



**Product Categories:** [Analytical Lab Instruments](#)

**Product Tags:** [Analytical Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/digital-potentiometer-2/>

### Product Description

DIGITAL POTENTIOMETER - Analytical Laboratory Equipment - Redox Potential. - Measurement display 4 digit LED range 0 to  $\pm 1999$  mV. - Accuracy  $\pm 0.1\%$  Power supply 230V $\pm$  10% 50Hz.

supplied with Glass, References (Calomel), Silver & Platinum Electrodes.

Range

pH

0 to 200

mV

0 to 2000

## Accuracy

pH

-----

mV

$\pm 1$  mV  $\pm 1$  digit

Resolution

0.1 mV

Slope correction

80% to 120%

Ambient Temperature

0° C to 50° C

Display

3½ digit LED display

Power

230 V  $\pm 10\%$  50Hz AC

Dimensions

275x175x85 mm

Accessories

- 1) Magnetic stirrer
- 2) Glass, Reference , Silver and Platinum Electrodes
- 3) Dust cover
- 4) Instruction manual.

Product data have been exported from - Labappara scientific instrument

Export date: Wed Nov 27 5:41:45 2024 / +0000 GMT