

DIGITAL POTENTIOMETER



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/digital-potentiometer-2/>

Product Description

DIGITAL POTENTIOMETER - Analytical Laboratory Equipment - Redox Potential. - Measurement display 4 digit LED range 0 to ± 1999 mV. - Accuracy $\pm 0.1\%$ Power supply 230V \pm 10% 50Hz.

supplied with Glass, References (Calomel), Silver & Platinum Electrodes.

Range

pH

0 to 200

mV

0 to 2000

Accuracy

pH

mV

± 1 mV ± 1 digit

Resolution

0.1 mV

Slope correction

80% to 120%

Ambient Temperature

0° C to 50° C

Display

3½ digit LED display

Power

230 V $\pm 10\%$ 50Hz AC

Dimensions

275x175x85 mm

Accessories

1) Magnetic stirrer 2) Glass, Reference , Silver and Platinum Electrodes

3) Dust cover

4) Instruction manual.

